


Diamond Ortho-strips

金剛石鋸弓



Creating additional space by interproximal enamel reduction (IPR) :

When opening the interdental space and reducing, finishing and polishing enamel in orthodontics, care must be taken to avoid grooves and scratches. It's important to proceed step by step, following a coarse-to-fine sequence.



Separating

OS40

0.40 mm



Reducing

OS30

0.30 mm




Contouring

OS25

0.25 mm

Item No.  1


Thickness (mm)



Finishing

OS20

0.20 mm



Pre-polishing

OS15

0.15 mm



Polishing

OS10

0.10 mm



No-Diamond Coated
with serrated edge

OS05

0.05 mm

